



HealthAtlas

Detailed User Manual Local Apple Health visualisation for macOS

English

Privacy-first local Apple Health visualisation for macOS

HealthAtlas - User Manual

Welcome

HealthAtlas reads an Apple Health export that you choose locally on your Mac. You then decide exactly which recognised data types appear in the dashboard. The app does not diagnose conditions or recommend treatment.

- Every normal launch starts without imported health data.
- Imported data stays in memory for the current session only. It is not uploaded or sent to an account.
- Language, theme, card count and data-type selection are stored separately for the Dev, Beta and Final app variants.
- While no personal data is loaded, the start screen also offers the matching GitHub manual and optional AI help. It copies only a general question with the public manual link to the clipboard; health values and local data stay in the app.
- Requirement: macOS 26 or later.


Important

HealthAtlas describes values and trends. It is not a substitute for professional medical advice, examination or diagnosis.

Quick start in 5 steps

The repository includes a fully synthetic demo export for a safe first test. It contains no personal health information.

- Open HealthAtlas. The empty overview presents a central import button.
- Click "Import ZIP or Export.xml ...".
- Choose Demo/AppleHealthDemo/Export.xml or your own Apple Health export.
- Open Sources after import and turn the data types you want on or off.
- Explore the selected data in Overview, Trends and Insights.



Empty start screen. The shown interface contains no health information.

Importing an Apple Health export

HealthAtlas accepts exactly two local formats: a direct Export.xml file or an Apple Health ZIP archive containing Export.xml. You do not need to unpack the ZIP first; HealthAtlas reads Export.xml directly from the archive.

- Files must be between 1 byte and 100 MB.
- Only one file can be selected at a time.
- The clinical companion file in an Apple Health export is intentionally not imported.
- There is no direct HealthKit connection and no cloud import.
- For a wrong ZIP, unreadable XML or an oversized file, HealthAtlas explains the issue and imports nothing.

Export on iPhone

In Apple Health, open Summary, tap your picture or initials, then select "Export All Health Data". Apple can change exact interface labels. Source: Apple Support, support.apple.com/en-in/guide/iphone/iph5ede58c3d/ios

Navigation and status


The frosted sidebar remains visible in every theme. It contains all five sections and shows "Private - Local only" at the bottom as a permanent privacy reminder.

- Overview: cards for selected data types.
- Trends: an interactive timeline for one selected type.
- Sources: selection of all data types recognised in the import.
- Insights: a descriptive summary for one data type.
- Design Studio: language and appearance.
- The GitHub and Discord icons above the privacy status open the project page and community in the default browser.

Overview

After a successful import, Overview shows only the types enabled in Sources. Every card has a type-specific accent colour and graphic treatment. The displayed numbers use the most recent available daily value.

- Choose 4, 8 or 12 to control how many cards appear per page.
- Previous and next controls appear when more selected types exist than fit on one page.
- Steps, energy, distance and flights climbed are shown as daily sums. Other numeric types are shown as daily averages.
- Name, value, unit and date come from the selected export and follow the app language formatting.




Overview with synthetic demo data. Colours aid recognition; they do not assess health.

Sources

Sources is the central selection area. It lists every recognised data type with an on/off switch, number of samples and a summary value. The table scrolls vertically when an export contains many types.

- Show all enables every recognised data type.
- Show none disables all data types.
- The switch in each row immediately determines whether a type is available in Overview, Trends and Insights.
- For a first personal import, the first four recognised types are selected by default. A compatible earlier selection is reused when its types occur in the new import.




Sources with synthetic demo data. This selection controls the content of all other areas.

Trends

This section visualises one active data type as a line. Select a period first, then choose one of the enabled types in the menu.

- 7D: last 7 days.
- 30D: last 30 days.
- 3M: last 90 days.
- 1Y: last 365 days.
- Click a point to reveal its date and formatted value above the chart. The selected point receives a visible marker.
- If fewer than two values exist in the selected period, the chart reports that there are not enough values.




Trend with synthetic demo data. Points are interactive and reveal their value after a click.

Insights

Insights condenses one enabled data type into a local descriptive card. Choose the type from the menu. The view shows the latest daily value, its date, the change from the previous value and a short 14-day line preview.

- A positive or negative change only describes the difference from the preceding available daily value.
- The screen does not supply normal ranges, alerts, ratings or medical conclusions.
- If only one dated value is available, the app states this clearly.




Insights with synthetic demo data. The card is a descriptive summary, not a diagnosis.

Design Studio

Use Design Studio to set the interface language and appearance. Changes take effect immediately and are stored locally for the current app variant.

- Language: Deutsch or English. Navigation, labels and known data-type names change with the selection.
- Clear Glass: a shared frosted layer with a calm cyan, blue, violet and pink glow plus subtle light points across the whole app.
- Midnight Glass: a dark blue glass surface.
- Aurora: a teal variation.
- Warmpaper: a warm red-violet variation.
- When macOS Reduce Motion is enabled and while importing, Clear Glass reduces or pauses motion.



Design Studio. Themes change appearance only, never the health data.

Privacy and limits

HealthAtlas is designed as a local visualisation. It has no account, analytics, advertising, tracking or hidden upload. The app does not send imported health values to HealthAtlas, GitHub, Discord or another service.

- Data remains in memory while the app is open and is not loaded again at the next normal launch.
- GitHub and Discord are opened only by an explicit click on the lower sidebar icons. They receive no imported values.
- ChatGPT, Gemini or Claude open only after an explicit click on the empty start screen. Before that, HealthAtlas only copies a fixed general introduction question with the public, language-specific manual link to the clipboard. Only Cmd+V lets the person decide whether to paste the question into that service.
- The bundled demo is synthetic and allows feature testing without personal data.
- HealthAtlas is not medical software and not emergency assistance.

Symptoms or uncertainty

Please contact qualified healthcare professionals. A visual change in a chart is not a medical statement.

Troubleshooting and Gatekeeper

HealthAtlas is ad-hoc signed because the project has no Apple Developer account. macOS Gatekeeper may therefore warn when opening a Dev, Beta or Final build for the first time.

- In Finder, Control-click the app and choose "Open". Confirm Open in the following dialog.
- If necessary, go to System Settings > Privacy & Security and choose "Open Anyway" for that exact HealthAtlas build.
- Do not disable Gatekeeper system-wide. Open only your own build or one from the official HealthAtlas project.
- For import errors, confirm that the ZIP contains Export.xml, the file is not empty, is below 100 MB and is stored locally.
- If no values appear, enable at least one type in Sources. Trends need at least two daily values in the chosen period.

Build variants

Dev, Beta and Final use separate local preferences. A theme or selection in Dev does not alter the settings of a Beta or Final app.

Complete control reference

This table summarises every current area, control and outcome.

Area	Control	Outcome
Header	Theme menu	Changes appearance immediately.
Header	Import	Opens ZIP or Export.xml picker.
Overview	4 / 8 / 12	Sets cards per page.
Overview	Arrows	Moves through additional selected cards.
Trends	7D / 30D / 3M / 1Y	Limits displayed days.
Trends	Data type menu	Changes the displayed timeline.
Trends	Data point	Shows point date and value.
Sources	Row switch	Enables or removes one type.
Sources	Show all / none	Enables or disables all types.
Insights	Data type menu	Changes the descriptive summary.
Design Studio	Language	Switches Deutsch and English.
Design Studio	Theme cards	Selects Clear Glass, Midnight Glass, Aurora or Warmpaper.
Lower sidebar	GitHub / Discord	Opens the external project or community link in the browser.